

S 3250 UV/VIS DETECTOR

The **Sykam S 3250 UV/Vis Detector** is a variable wavelength UV/Vis detector for routine analysis and sophisticated research. The dual lamp design offers a wavelength range of 190 – 800 nm with a low baseline noise. The front-accessible flowcell can easily be exchanged, as can be the lamps which are accessible through a side panel in the instrument housing.



■ Integrated Wavelength Program

The **S 3250 UV/Vis Detector** features a wavelength program to change the selected wavelength over time. With this feature the optimum wavelength can be selected for each analyzed substance according to its retention time.

■ Integrated Peak Detector

The integrated Peak Detector works as a basic fraction collector. The peak detection level can be freely programmed for peak start and peak end to enhance the collection purity. An integrated 24V output for switching a solenoid valve is used for the fraction collection, which is automatically operated with a selectable time delay.

■ Optional – Dual-Wavelength

The **S 3250 UV/Vis Detector** is available with an optional second wavelength. This feature enhances the Wavelength Program feature that you can measure 2 different wavelengths at the same time. A second D/A converter output comes with this option to keep the system flexible to be used with any data acquisition software available.

■ Optional – Online-Scan

Another option for the **S 3250 UV/Vis Detector** is the Online Scan. With the Online Scan whole spectrum information can be gathered at a certain time. This scan information is stored internally and can be accessed at any time. The Online Scan is a good alternative to a full UV PDA detector.




Figure: S 3250 Peak Detector




Figure: S 3250 Optical Module

■ Technical Specifications*

| | |
|------------------------------------|--|
| Wetted Materials: | Stainless Steel / PEEK* |
| Baseline Noise: | $\pm 1 \times 10^{-5}$ AU (@240 nm, 1 sec. Risetime) |
| Baseline Drift: | $<3 \times 10^{-4}$ AU/h |
| Wavelength Range: | 190 – 800 nm |
| Wavelength Accuracy: | ± 2 nm |
| Linearity: | > 2.0 AU |
| Light Source | Deuterium Lamp, Tungsten Lamp |
| Wavelength Program: | Programmable, 10 steps |
| Analog Output: | 1x 1 V (optional: 2x 1V) |
| Control Features: | Internal Peak Detector with +24 V solenoid switching output. |
| Dimensions: (W x H x D) | 310 x 165 x 478 mm |
| Power Supply: | 100 - 250 ~V (47 - 63 Hz) |

* depending on configuration

■ Order Information

| Catalog No | Description |
|------------|---|
| 10 31 008 | S 3250 UV/Vis Variable Wavelength Detector |
| 10 31 009 | option: Dual Wavelength |
| 10 31 010 | option: Online Scan |
| 10 32 009 | S 3250 Flowcell, analytical, Stainless Steel |
| 10 32 010 | S 3250 Flowcell, analytical, PEEK |
| 10 32 011 | S 3250 Flowcell, micro, Stainless Steel |
| 10 32 012 | S 32505 Flowcell, micro, PEEK |
| 10 32 013 | S 3250 Flowcell, preparative, Stainless Steel |
| 10 32 014 | S 32505 Flowcell, preparative, PEEK |
| 40 10 002 | Spare Deuterium Lamp |
| 40 10 001 | Spare Longlife Tungsten Lamp |

Sykam GmbH

Systeme & Komponenten analytischer Meßtechnik

Gewerbering 15
86922 Eresing
Germany

Tel.: ++49 (8193) 93 82 - 0
FAX: ++49 (8193) 93 82 - 20
EMail: Info@sykam.com
Web: http://www.sykam.com

* All technical specifications may be subject to change.